

# IMSPECTOR

## imaging spectrographs

### General Description

Middleton Research offers the complete line of Specim Ltd. ImSpector™ Series imaging spectrographs for spectral imaging applications in both industry and scientific research for every wavelength region. Spectrographs are complex optical components that change a camera to a hyperspectral imaging device with full contiguous spectral information and high spectral resolution. The spectrograph can be attached and integrated with a monochrome camera using a standard C- or U-mount to build a line imaging hyperspectral camera. The spectrograph may be integrated by the end-user. Specim Ltd. also offers an integration service; contact Middleton Research for more information.

The unique design of the transmission imaging technology makes the spectrograph an inexpensive, light-weight, rugged component designed to withstand industrial use and harsh environments. The ImSpector Series spectrographs combine high diffraction efficiency over a broad spectral range. The transmission grating-based optical design allows excellent image quality and a short focal length. The spectrograph designs are also independent of the polarization state of the incoming light. The spectrographs can be integrated with multi-channel fiber optic arrays, creating devices capable of simultaneous spectral measurement of up to 120 points.

### Accessories


Various accessories can be combined or integrated with the spectrographs, including fore optics (special hyperspectral lenses), fiber optics, or filters to improve the quality of measured data. Hyperspectral fore optics are available for wavelength regions between 400 nm and 12,000 nm. Multiple point fiber optics are offered for wavelength regions between 200 nm and 2,500 nm. Find more details and options in the Accessories chapter.

### Slits

Each ImSpector model has a standard input slit with a fixed length. The length of the slit defines the length of the usable spatial axis (i.e. sensor width) and thus also defines the spatial resolution. Together with the chosen lens, the length of the slit also determines the field-of-view.

The slit length is fixed in a spectrograph, however a range of standard slit widths are available for each ImSpector and must be defined at the time the spectrograph is ordered so that the installation and calibration can be performed at the factory. All default and standard slit widths are listed in the specifications tables for each ImSpector. Different slit widths can be manufactured, so a non-standard slit width may be chosen for an additional fee.

The slit width affects three parameters: spectral resolution, throughput, and the width of the imaged line. The spectral resolution can be improved by choosing a narrower slit width. However, please note that the pixel size of the focal plane detector and other optical design factors also affect the spectral resolution. The optical throughput of the spectrograph is also affected by the width of the slit. Decreasing the slit width will cause reduced throughput, which can be offset by an increase in the integration time or a brighter illumination source. The width of the slit (together with the lens) determines the width of the imaged line.



*MWIR M50M ImSpector shown with lens*

## Selecting a Spectrograph

When choosing a spectrograph, the most important factors to consider are wavelength range and resolution. Minimizing aberrations is essential for high resolution applications. For more information on keystone and smile aberrations, please refer to the glossary at the end of the catalog.

The ImSpector spectrograph product line includes four different series: Enhanced, Standard, Fast, and M-Series.

- The Enhanced E-Series offers the best optical performance, and is recommended for high resolution applications (< 3 nm). Models are available for UV, visible, Raman, VNIR, NIR and SWIR ranges and are suitable for large area sensors. The enhanced series offers higher light throughput than the Standard Series and minimizes keystone and smile aberrations to sub-pixel levels.
- The Standard Series contains basic spectrographs used mainly for machine vision and other, less spectroscopically-demanding applications. Models are available for VIS and VNIR and for 2/3" or 1/2" sensors. The spectral resolution ranges from 6 - 11 nm.
- The Fast Series spectrograph for the VNIR region is designed for both high speed (> 200 fps) and low light level applications. The fast spectrograph is used to achieve short read-out time (high frame rate) and high light throughput (short integration time) at a lower spectral resolution (15 nm).
- The M-Series is a new line of spectrographs for the VNIR, MWIR, and LWIR ranges that produces the best image quality by providing an increase in spatial resolution. These spectrographs use a transmissive dispersive element based on a new optical design. As a result, the M-Series spectrographs are suitable for use with larger detectors while remaining lightweight and compact.

Wavelength Region	Spectral Range	Model	Spectral Resolution
Raman	530 – 630 nm	R6E	0.3 nm
	770 – 980 nm	R10E	0.3 nm
UV	200 – 400 nm	UV4E	2 nm
VIS	380 – 800 nm	V8H	8 nm
		V8	6 nm
		V8E	2 nm
VNIR	400 – 1000 nm	V10H	11.2 nm
	400 – 1000 nm	V10	9 nm
	400 – 1000 nm	V10E	2.8 nm
	400 – 1000 nm	Fast10	15 nm
	350 – 1000 nm	V10M	1.5 nm
NIR	900 – 1700 nm	N17E	5 nm
SWIR	1000 – 2500 nm	N25E	8 nm
MWIR	3000 – 5000 nm	M50M	35 nm
LWIR	8 – 12 $\mu$ m	*	*

\* LWIR spectrographs are available only as full hyperspectral imaging cameras. Refer to Cameras chapter for details.